



accutest® Uriscreen™

Rapid Screen For Urinary Tract Infection (UTI)

Product Overview



The Accutest® Uriscreen™ is a rapid screening test which detects both bacteriuria and pyuria (the presence of somatic cells in urine) by identifying catalase activity (catalase is an enzyme present in many eukaryotic cells).

This test provides **physicians' offices, clinical labs, long-term care, urgent care and other outpatient settings** with a rapid and accurate screening method for urinary tract infections.

Accutest® Uriscreen's **negative predictive value (NPV) of 95% is higher than the NPV of the urine dipstick.** This test is a cost-effective screening procedure which can be used to complement conventional culture methods.

Ordering Information

Catalog #.	Product Name	CPT Code
ID601	Accutest® Uriscreen™	81007QW
ID104-01	Accutest® Uriscreen™ Controls	



accutest[®] Uriscreen[™]


Rapid Screen For Urinary Tract Infection (UTI)

Highlights & Benefits

- ✓ **ECONOMICAL**
Cost effective; can be used to complement conventional culture methods.
- ✓ **EASY TO USE**
Rapid results in 2 minutes or less; requires no additional instrumentation.
- ✓ **LONG SHELF LIFE**
Thirty six (36) months at room temperature from date of manufacture.
- ✓ **PERFORMANCE**
Highly sensitive test for bacteriuria and the presence of somatic cells in urine.
95% Negative Predictive Value (higher than urine dipstick method)

REIMBURSEMENT
CPT CODE: 81007QW


Simple Procedure



1.5 – 2 mL of urine specimen is added to a test tube



4 drops of hydrogen peroxide solution are added to the tube and mixed



The appearance of white foam indicates a positive result



NEGATIVE



POSITIVE

Indicates presence of UTI (bacteriuria/hematuria/pyuria)

Kit Contents

20 TESTS PER KIT

- 1 dropper bottle with 10% hydrogen peroxide (H₂O₂) solution
- 20 stoppered tubes with test reagent powder
- 20 disposable 20mL pipettes
- 1 instruction manual

For more information visit our website at
www.JantDx.com/Accutest-Uriscreen



JANT PHARMACAL CORPORATION
LABORATORY & POINT OF CARE DIAGNOSTICS
www.JantDx.com

JANT PHARMACAL CORPORATION
16530 VENTURA BLVD. STE 512 ENCINO, CA 91436
PHONE: 800.676.5565 FAX: 818.986.0235
EMAIL: info@jantdx.com