



Status™

lifeSign

COVID-19/Flu A&B

A Rapid Immunoassay for the Simultaneous Direct Detection and Differential Diagnosis of SARS-CoV-2, Influenza Type A and Type B Antigen from anterior nasal and nasopharyngeal swab specimens

Infections with these viruses may present similar symptoms. Can you tell them apart? **We Can.**




1 test + 1 specimen = 3 Results




- **COVID-19 - Anterior Nasal - Sensitivity 93.8%, Specificity 100%**
Nasopharyngeal - Sensitivity 93.1%, Specificity 100%
- **Flu A - Sensitivity 91.4%, Specificity 95.7%**
- **Flu B - Sensitivity 87.6%, Specificity 95.9%**
- **FDA Emergency Use Authorization (EUA)**
- **Visually read in 15 minutes**
- **Flocked swabs for superior specimen collection and patient comfort**

ORDERING INFORMATION

PRODUCT	# OF TESTS	ITEM #
Status™ Covid-19/Flu A&B	25	ID33225



Easy 3 Step Procedure

1		Tear the tab off the Extraction Reagent Capsule and dispense entire contents into the Extraction Well.
2		Insert the specimen swab in the Swab Stand. • Rotate swab 3 times to mix the specimen. • Let stand 1 minute. • Rotate swab 3 times again and discard the swab.
3		Raise the device upright and let stand 1–2 seconds. Gently tap the device to ensure the liquid flows into the hole. Lay the device back down.
Set a timer for 15 minutes.		
4	Read test results at 15 minutes. NOTE: False positive or false negative results can occur if the test is not read between 15 and 20 minutes.	

- Refer to package insert and other training materials for complete test procedure details.
- Status™ COVID-19/FLU A&B is an antigen test for use with anterior nasal and nasopharyngeal swab specimens.
- Status™ COVID-19/FLU A&B test has been granted emergency use authorization (EUA) by the FDA for use at the point-of-care setting. During the duration of the emergency this test can be performed in a patient care setting that is operating under a CLIA Waiver, Certificate of Compliance or Certificate of Accreditation.

 Proudly made in the USA