Sweat Control for Cystic Fibrosis Testing / Level 1, Level 2, Level 3

LOT 229141
229142
229143

REF 1190-01

2028-05-31

European Conformity
CE-Konformitätsbescheinigung
Conformité aux normes européennes
Conformità europea
Conformidad europea

Lot Number
Designation du lot
Numero di lotto
Denominación de lote

Manufactured by
Hergestellt von
Fabricato per
Fabricado por

In vitro diagnostic use
In-vitro Diagnosticum
Pour diagnostic in vitro
Per uso diagnostico in vitro
De uso diagnóstico in vitro

Biological Risk
Biosicherheit
Risque biologique
Rischio biologico
Peligro biológico

Contents of kit
Inhalt der Packung
Contenu du coffret
Contenido della confezione
Contenido del estuche

Catalog No.
Bestellnr.
N° de catalogue
Catalogo n.
N° de catálogo

Caution, See Product Insert
Achtung, Siehe Packungsbeilage
Attention, voir notice d'utilisation
Attenzione, vedere il foglietto illustrativo del prodotto
Atención, consulte el folleto del producto

Authorized Representative
Bevollmächtigter
Représentant agréé
Rappresentante autorizzato
Representante autorizado

Temperature limitation
Temperaturbegrenzung
Limites de température
limite di temperatura
límite de temperatura

Consult instructions for use
Gebrauchsanweisung beachten
Consulter les instructions d'utilisation
Consultare le istruzioni l'uso
Consulte las instrucciones de uso

Use by (last day of month)
Verwendbar bis letzter Tag des Monats
Utilisable jusqu'à (dernier jour du mois indiqué)
Da utilizzare prima del (ultimo giorno del mese)
Estable hasta (último día del mes)

English

Intended Use

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are assayed controls intended as a means of monitoring the measurement of conductivity and osmolality by cystic fibrosis screening procedures and to validate the analysis of Chloride, Potassium, and Sodium by confirmatory methods used to quantitate patient samples.¹

Product Description

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing are liquid, ready to use, requiring no reconstitution or dilution. They are supplied in three levels, nine vials per box, 3 x 3 mL each level. They are prepared in an aqueous simulated human sweat matrix² to which preservatives have been added.

Warnings and Precautions

The Quantimetrix Sweat Controls for Cystic Fibrosis Testing do not contain human source material. No special precautions are required in handling the product other than those routinely used in the laboratory.

Warning ⚠ Hazard (H) and Precautionary (P) Statements

Contains Mixture, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone. Urea, Sodium Chloride, Potassium Chloride.

H317 – May cause an allergic skin reaction

P261 – Avoid breathing vapors, mist, or spray.

P272 – Contaminated work clothing should not be allowed out of the workplace.

P280 – Wear protective gloves, protective clothing, and eye protection.

P302+P352 – IF ON SKIN: Wash with plenty of water.

P333+P313 – If skin irritation or rash occurs: Get medical advice/attention.

P362+P364 – Take off contaminated clothing and wash it before reuse.

P501 – Dispose of contents/container in accordance with local, regional, national, and international regulations.

Safety Data Sheet (SDS) available for professional users at quantimetrix.com

Storage and Stability

The Sweat Controls should be stored refrigerated (2°C– 8°C) or at room temperature (18°C– 25°C). Do not freeze. When stored unopened at 2°C– 8°C or at 18°C– 25°C the controls are stable until the expiration date stated on the label. When stored at 2°C– 8°C between each use the controls are stable for 3 months after opening. Discard the controls if turbid or any evidence of microbial contamination is present.

Procedure

Remove the controls from the refrigerator and allow them to come to room temperature (18°C– 25°C), at least 15 minutes, depending on remaining volume. Invert gently to ensure homogeneity of the contents. Immediately recap the control when not in use.

Expected Values

The expected values have been established in the Quantimetrix laboratory and from interlaboratory data using the listed manufacturer's instruments. Use of control materials having known component concentrations are an integral part of diagnostic procedures. Monitoring of control values establishes intralaboratory parameters for accuracy and precision of the test method. Each laboratory should establish its own quality control standards for the test method used. A quality control program, Quantrol, is available from Quantimetrix. Please inquire.

Limitations

Any future changes made by the manufacturer of a test method may give a different value from the one indicated. Limitations of the test method are included in the information provided by the instrument manufacturer. Technical updates can be found on our website.

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Gemisch, 3(2H)-isothiazolone, 5-chloro-2-methyl- with 2-methyl-3(2H)-isothiazolone. Urea, Sodium Chloride, Potassium Chloride.

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Quantimetrix Corporation
MDSS GmbH