



# NicAlert

## **Urine/Saliva Cotinine Test** **For Accurate Measurement of Nicotine** **Consumption & Tobacco Exposure**



- Tests for cotinine (nicotine metabolite) in urine and saliva
- One-step test with results in minutes
- Semi-quantitative test for accurate measurement of nicotine and tobacco consumption
- Detects first and second-hand smoke
- Identifies 6 levels of tobacco exposure
- Highly accurate point of care test to determine level of exposure to tobacco products

**DETECTS  
SECOND-  
HAND  
SMOKE\***



### **Advantages:**

- NO mixing or measuring of reagents or calibration of equipment required
- NO instrumentation or trained personnel needed to read results
- NO "first morning" urine required
- NO cross-reactivity with drugs of abuse and other compounds

\* With urine test only




accutest®

# NicAlert

## DETECTS SECOND-HAND SMOKE

### Easy Procedure:

1. Dip test strip as instructed.
2. Place on surface below.
3. Wait for blue line to disappear.
4. Number of lowest red band is result.



CAT. # DS43: NICALERT SALIVA  
CAT. # DS47NY150: NICALERT URINE

Ideal for use in:

### Research and Pharmaceutical Clinical Trials:

- Identifies smoking status
- Smoking can induce liver enzymes, influence drug pharmacokinetics and rate of drug elimination from body

### Smoking cessation programs:

- Readily demonstrates biological effect to the smoker
- Shows the effect on others via second-hand smoke
- Provides immediate feedback to doctor and patient
- Provides a measure for smokers trying to quit
- Provides a tool to initiate patient counseling

### Insurance risk assessments:

- Identifies levels of smoking/tobacco use
- Acts as an aid to individual risk assessment

### Importance for Healthcare Professionals:

- Identifies which patients are smoking
- Detects second-hand smoke exposure
- Measures level of nicotine dependency in smokers
- Assists patients in quitting
- Identifies smokers who relapse

### Important Patient Groups:

- Pregnant Women
- Smokers
- Asthmatics and their families
- Cardiac and Respiratory patients
- Teenagers



accutest®

DISTRIBUTED BY:

### JANT PHARMACAL CORPORATION

16530 Ventura Blvd., Suite 512, Encino, CA 91436  
800.676.5565 818.986.8530 Fax 818.986.0235  
www.accutest.net info@accutest.net

Accutest is a registered trademark of Jant Pharmacal Corporation