accutest H. pylori Urease Test

The Next Generation in Helicobacter Pylori Testing



- RAPID RESULTS
- NO REFRIGERATION
- NO REAGENTS



RAPID UREASE TEST FEATURES:

FEATURE	BENEFIT
No Refrigeration Required	Can be stored at room temperature No need to keep temperature logs.
Rapid Results	Positive in as little as 1 minute, negative results in 1 hour
Accurate	99% Sensitivity and 97% Specificity as shown in major US studies comparing Accutest Rapid Urease Test to histology and other R.U.T.'s
Dry Test	No incubation or reagents required
Cost Effective	Low cost per test
Simple and Convenient	CLIA waived
2 Year Shelf Life	Good stability, superior shelflife

5 SIMPLE STEPS

1. Peel back label


2. Place biopsy on test


3. Re-seal test


4. Wait up to one hour


5. Read results













JANT PHARMACAL CORPORATION


16530 Ventura Blvd., #512, Encino, CA 91436 **800.676.5565 818.986.8530** Fax 818.986.0235 www.accutest.net info@accutest.net

accutest H. pylori Urease Test

The Next Generation in Helicobacter Pylori Testing



- RAPID RESULTS
- NO REFRIGERATION
- **NO REAGENTS**





- SIMPLE
- COST EFFECTIVE
- RESULTS IN 1 **HOUR OR LESS**

COMPARISON TO OTHER DIAGNOSTIC METHODS:

- Much more expensive
 - Not timely

HISTOLOGY:

- Accuracy depends on experience and
- appropriate stains Relatively expensive

- Does not define current H. pylori infection Less accurate than histology or Accutest
- H. pylori Urease Test

Product # Description Packaging

ID750 Accutest H. pylori Urease Test 50 Each

CPT Code: 87077 QW

IDC73 Urease/Glycerin Control Solution 3mL Dropper Bottle

Positive Control for Accutest h. Pylori Urease Test



JANT PHARMACAL CORPORATION

16530 Ventura Blvd., #512, Encino, CA 91436 800.676.5565 818.986.8530 Fax 818.986.0235 www.accutest.net info@accutest.net