

YOUR RESOURCE FOR POINT OF CARE TESTING


Jant Pharmacal's extensive menu of CLIA-Waived tests helps you
MAKE AN IMPACT IN PATIENT CARE!



accutest[®]

Point of Care Diagnostics

CLIA-Waived



JANT

PHARMACAL CORPORATION
Laboratory & Point of Care Diagnostics
WWW.JANTDX.COM

INFECTIOUS DISEASE



AccuStrip® Value + Rapid Strep A Strip Test (Catalog # ID513, 25 test strips/box)

- Results in 5 minutes
- High Sensitivity and Specificity
- CLIA-Waived
- Room temperature storage



Accutest® Value+ Mononucleosis Rapid Test (Catalog # ID516, 25 tests/box)

- FDA cleared with no age restrictions
- Results in 8 minutes
- 99.9% Sensitivity and 98.8% Specificity
- Room temperature storage
- CLIA-Waived for whole blood



Accutest® H Pylori Rapid Test (Catalog # ID392, 20 tests/box)

- Results in 3-7 minutes
- Room temperature storage
- CLIA-Waived for whole blood



THYROID

Accutest® TSH Test for Hypothyroidism (Catalog # MD801, 20 tests/box)

- Requires only one drop of blood from finger-stick
- Results in 10 minutes
- 98.5% Sensitivity, 96.9% Specificity
- CLIA-Waived
- Room temperature storage



PREGNANCY

Accutest® Value+ Urine hCG Tests

- Available as test strips (Catalog # PF851, 25 tests/box) or cassettes (Catalog # PF853, 25 tests/box)
- Results in 5 minutes or less
- >99% Sensitivity
- Cut-off: 25mLU.mL
- CLIA-Waived

Accutest® Value+ Serum/Urine hCG Tests

(Catalog # PF864, 25 tests/box)





DRUGS OF ABUSE

Accutest® Drug Test Cups

- Self-contained
- Rapid Results in 5 minutes
- Leak-proof cup may be sent to lab for confirmation
- Tests up to 14 different drugs (Amphetamine, Barbiturates, Benzodiazepines, Buprenorphine, Cocaine, Ecstasy, Marijuana, Methadone, MDMA, Methamphetamine, Opiates, Oxycodone, Phencyclidine, Tricyclic Antidepressants) and 6 adulterants
- Tests available as dips or cassettes
- CLIA-Waived
- Room temperature storage



CANCER SCREENING

Accutest® Single or Dual iFOB Tests

- Sensitivity 99.2%, Specificity 96.7%, Accuracy 98.3%
- No dietary or drug restrictions
- Easy to Use collection materials
- Extended Stability of sample in buffer solution
- Cutoff is 50ng hHb/ml

URINALYSIS

AccuStrip® URS Reader and 10 Parameter Urine Reagent Strips (Catalog # UA850, URS Reader)

- Touch Screen Operation
- Battery and External Power Options
- Memory for 200 Measurements
- Reads 50 strips/hour
- Accustrip® URS 10 parameter strips test for Blood, Urobilinogen, Bilirubin, Protein, Nitrite, Ketone, Glucose, pH, Specific Gravity and Leukocytes (Catalog # UA870, 100 strips/box) and are CLIA-Waived for visual or analyzer testing
- Promotional Reader & Strip Package (Catalog # PRM804, 1 URS reader and 4 bottles of strips)



AccuTest® URS 10 Parameter Urine Reagent Strips (Catalog # UA710A, 100 tests/bottle)

- CLIA-Waived when read visually
- Reliable Sensitivity
- Easy to Read



Accutest® UriScreen™ (Catalog # ID601, 20 tests/box)

- Rapid screen for urinary tract infections
- Results in 2 minutes
- High Sensitivity, 95% Negative Predictive Value
- CLIA-Waived



CHOLESTEROL AND GLUCOSE

LipidPlus® Total Cholesterol and Glucose Testing System with portable hand-held battery operated analyzer and printer.

(Catalog # MD532, 10 Total Cholesterol tests/box)

(Catalog # MD540, 50 Glucose strips/box)

- Results in less than 2 minutes
- Small Sample Size: 5µl capillary whole blood for Total Cholesterol and 0.3 µl for Glucose
- CRMLN Certified for Total Cholesterol
- CLIA-Waived
- Room temperature storage



COAGULATION

Accutest® Coag-Sense™ PT/INR Monitoring System

(Catalog # MD03P60, PT/INR System)

(Catalog # MD03P56, 25 tests/box)

- Direct Clot Detection Technology
- Accurate results in less than one minute
- Small Sample Size: 10µl capillary whole blood
- Room temperature storage
- Portable analyzer included with purchase of 4 boxes of professional strips (PRM PT/INR)
- CLIA-Waived



SPOTCHEM™ EZ CLIA-WAIVED CHEMISTRY ANALYZER

This small benchtop chemistry analyzer uses dry chemistry technology. Testing is performed on reagent strips with results available within 3-15 minutes. The menu consists of 21 tests with a throughput of 48 tests per hour. The instrument accepts whole blood and serum samples and is calibrated automatically with the swipe of a magnetic card. Simply add sample and insert centrifuge cup, place test strip and press start. At the conclusion of testing, results are reported via the SPOTCHEM EZ printer or can be downloaded via LIS system.

(Catalog # SPEZ443, Analyzer)

CLIA-Waived tests include the following:

Individual Tests: Glucose, ALT, AST, Total Cholesterol, Triglycerides, HDL, BUN, Creatinine, Total Protein

Lipid Panel: Total Cholesterol, Triglycerides, HDL, Calculated LDL

General Chemistry Panel 1: ALB, BUN, GLU, CRE, CAL

General Chemistry Panel 2: ALP, T-BIL, ALT, AST, TP

