



## SAFETY INFORMATION

ON

**Cat#: PF418**  
**ACCUTEST® pH**  
**PHENAPHTHAZINE PAPER**

Accutest pH test paper is classified as an article under OSHA regulations. Accutest pH test paper does not require formal material safety data sheets or certificates of analysis (under osha standard – “hazardous communications” – 29 cfr 1910.1200). It is safe for shipment by air, sea or ground transportation.

The test paper consists of indicator dyes that have been impregnated into the paper. The paper is to be used for in vitro diagnostic use only and not to be used internally. The concentration of impregnated dye is small.

Accutest pH test paper consists of stable, non-reactive materials. It is non-volatile and not radioactive. Accutest pH test paper and its components do not contain latex, mercury, DEHP or lead. The paper can be discarded according to local environmental waste procedures and regulations.

Accutest pH test paper has a shelf life of 3 years and should be stored in a dry protected area away from sunlight.

For any additional information, please contact:

[www.JantDX.com](http://www.JantDX.com)  
**Jant Pharmacal Corporation.**  
Encino, California 91436, USA  
Tel: 800-676-5565; Fax: 818-986-0235  
info@jantdx.com