

## ACCUTEST Drug Test Cup 5 with Adulteration Screen

- Hygienic self-contained cup format.
- Easy collection.
- Easy to read results.
- No reagents needed.
- Stores at room temperature.
- Built-in controls to validate results.
- Adheres to SAMHSA cutoff levels
- Self activated test - Just add urine and read results.
- **Tests for:**
  - Amphetamines
  - Cocaine
  - THC (marijuana)
  - Opiates
  - PCP



Adulteration Screen tests for: Creatinine, pH, Specific Gravity, NBP

- **Catalog # DS854625**
- 25 tests/box