

JANT PHARMACAL CORPORATION


LABORATORY & POINT OF CARE DIAGNOSTICS

ValuPak™ Strep A Tests

Accustrip® ValuPak™ Strep A Test Strip is a one-step immunoassay test for the rapid, qualitative detection of Group A Streptococcal antigen from throat swabs.



- ✓ Easy-to-use & CLIA Waived
- ✓ High sensitivity and specificity
- ✓ Rapid results within 10 minutes



50 TESTS PER KIT

- CAT # ID350
- 50 Test Strips
- 50 Extraction Tubes
- 50 Swabs
- CPT Code: 87880QW

ValuPak™ Series- Quality and Convenience for Less!



JANT PHARMACAL CORPORATION

16530 VENTURA BOULEVARD, SUITE 512, ENCINO, CA 91436
PHONE: 800.676.5565 818.986.8530 FAX: 818.986.0235

info@jantdx.com www.jantdx.com