



# New controls for serology

Get more done with consolidated QC solutions for infectious serology

Thermo Scientific™ MAS™ Omni•Infectious™ Quality Controls for viral blood stream infections allow you to optimize your daily workflow and increase your productivity. The true **third-party positive and negative control panels** for hepatitis, HIV, and HTLV are consolidated, allowing you to save time. Furthermore, you can assess the performance of your serological assays in an objective and unbiased manner as opposed to working with in-kit controls.

**#MASproductivity**  
continues with serological controls



**30 DAYS**  
open vial stability @ 2–8 °C\*



**READY TO USE**  
straight from the fridge



**HUMAN**  
serum matrix



BSI  
Positive



Infectious  
Negative



Individual  
Positive Controls

### Consolidated controls for efficiency






- ✓ **HIV** antibodies to human immunodeficiency virus 1/2
- ✓ **HTLV** antibodies to human T-lymphotropic virus I/II
- ✓ **HBs** Hepatitis B Surface antigen (nonreactive in negative control for HBsAg and anti-HBs)
- ✓ **HCV** antibodies to Hepatitis C Virus

### Independent use

- ✓ **HIV-p24** HIV p24 antigen
- ✓ **HBs** antibodies to HBsAg
- ✓ **TP (Syphilis)** antibodies to *Treponema pallidum*

\*Not all controls have 30 days open vial stability.

## Thermo Scientific MAS Quality Controls for infectious serology

Part Number	Product	Analyte	Pack Size	Shelf Life @ 2 – 8 °C	Open Vial Stability @ 2 to 8 °C	Matrix
 10027480	MAS Omni•Infectious BSI Positive Control Panel	anti-HIV 1/2 anti-HTLV I/II HBsAg anti-HBc anti-HCV	5 x 5 mL	13 months	60 days	Human serum
 10027484	MAS Infectious Negative Control	anti-HIV 1/2 anti-HTLV I/II HBsAg anti-HBs anti-HBc anti-HCV anti-TP	10 x 5 mL	24 months	60 days	
 10027481	MAS HIV p24 Positive Control	HIV p24	5 x 5 mL	24 months	30 days	
 10027482	MAS anti-HBs Positive Control	anti-HBs	5 x 1 mL	24 months	30 days	
 10027483	MAS Syphilis Positive Control	anti-TP (Syphilis)	5 x 1 mL	18 months	30 days	



MAS Omni•Infectious BSI Positive Control Panel



MAS Infectious Negative Control



MAS HIV p24 Positive Control



MAS anti-HBs Positive Control



MAS Syphilis Positive Control

Call your local Thermo Fisher Scientific representative or Tech Services for additional information.

### Clinical Diagnostics

**USA**  
46500 Kato Road  
Fremont, CA 94538 USA  
Tel. 800-232-3342 / +1 510-979-5000  
sales.diagnostics.fmt@thermofisher.com

**Asia Pacific**  
Tel: +65 6499 9999  
cdd.asia.info2@thermofisher.com

**Australia**  
Tel. +61 1800 333 110  
auinfo@thermofisher.com

**Canada**  
Tel. +1 905-286-4290  
CDD.Canada@thermofisher.com

**China**  
Tel. +86 400 650 5118  
sales.china@thermofisher.com

**France, Belgium (Brussels & Wallonia) & Luxembourg**  
Tel. +33 (0) 1 40 86 65 00  
cdx.fr.info@thermofisher.com

**Germany**  
Tel. +49 (0) 3302 883 608  
cdx.de.info@thermofisher.com

**Hong Kong**  
Tel. +852 3107 7600  
SDG.HK.Sales@thermofisher.com

**Italy**  
distributor.cdd@thermofisher.com

**Japan**  
Tel. +81 (0) 120 489 211  
JPYOK-CDD.QC@thermofisher.com

**Netherlands & Belgium (Flemish)**  
Tel. +49 3302 883 660  
Info.benlux.cdd@thermofisher.com

**New Zealand**  
Tel. +61 0800933 966  
nzinfo@thermofisher.co.nz

**Nordic**  
Tel. +46 8 594 691 30  
info.nordic.cdd@thermofisher.com

**Russia**  
Tel. +7 495 739 76 41  
info.btd.moscow@thermofisher.com

**South Africa**  
Tel. +27 11 792 6790  
support-za.idd@thermofisher.com

**Spain & Portugal (Iberia)**  
Tel. +34 93589 8338  
CDD.info.ES@thermofisher.com

**Switzerland & Austria**  
Tel. +41 26 663 86 70  
cdx.ch.info@thermofisher.com

**United Kingdom & Ireland**  
Tel. +44 1442 868 940  
cdx.UK.info@thermofisher.com

**For Countries Not Listed**  
distributor.cdd@thermofisher.com