




BENCHTOP TOXICOLOGY ANALYZER

IMMTOX™ 270

Small, simple
and smart

The ImmTox™ 270 benchtop analyzer offers comprehensive toxicology screening solutions for physician offices, treatment centers and independent laboratories.

- Compact design with easy-to-use interface
- Comprehensive test menu CLIA¹ categorized as moderate complexity
- Low cutoff concentrations
- Fast results
- High throughput
- Walkaway operating functionality
- Minimal maintenance



ImmTox™ 270

FEATURES AND BENEFITS

- 270 dual reagent tests per hour
- Reagent capacity: 24 assays
- Onboard patient sample capacity: up to 30
- Ideally sized for 200-1200 patient samples per month
- Intuitive, touch-screen user interface
- Bi-directional LIS communication
- Fully traceable results and quality control data
- Continuous access to samples without interruption of the testing process
- Reusable cuvettes and wedges
- Fast start-up
- Strong technical support

LIS COMMUNICATION SOFTWARE AVAILABLE

- Easy management of patient records
- Flexibility to generate summary reports for patients and quality controls
- EMR² integration options
- HL7³ interface

SPECIFICATIONS

- Dimensions: 32”(W) x 27”(D) x 21”(H) – 38”(H) with open cover
- Electrical: Standard 110V outlet
- Self-contained analyzer
 - Water consumption: 3.8 liters
 - External water or drainage connections not required
 - Easy to install and operate

TEST MENU

URINE DRUG SCREENING REAGENTS	CUTOFF LEVEL (ng/mL)	MODERATE COMPLEXITY
6-acetylmorphine (6-AM Heroin metabolite)	10	<input checked="" type="checkbox"/>
Amphetamine	500 or 1,000	<input checked="" type="checkbox"/>
Barbiturates	200	<input checked="" type="checkbox"/>
Benzodiazepines	200	<input checked="" type="checkbox"/>
Benzoylcegonine (Cocaine metabolite)	150 or 300	<input checked="" type="checkbox"/>
Buprenorphine	5	<input checked="" type="checkbox"/>
Cannabinoids (THC)	50	<input checked="" type="checkbox"/>
Carisoprodol	100	
EDDP (Methadone metabolite)	100 or 300	<input checked="" type="checkbox"/>
Ethyl Alcohol	100 mg/dL	
Ethyl Glucuronide (EtG)	500 or 1,000	
Fentanyl (SEFRIA™)	1	<input checked="" type="checkbox"/>
Fentanyl (HEIA™)	2	
Ketamine	100	
Meperidine	200	
Methamphetamine	500 or 1,000	<input checked="" type="checkbox"/>
Opiates	300 or 2,000	<input checked="" type="checkbox"/>
Oxycodone	100 or 300	<input checked="" type="checkbox"/>
Phencyclidine (PCP)	25	<input checked="" type="checkbox"/>
Synthetic Cannabinoids-1 (JWH-018, JWH-073, AM-2201 and its metabolites)	10	
Synthetic Cannabinoids-2 (UR-144, XLR11 and its metabolites)	10	
Synthetic Cannabinoids-3 (AB-PINACA, ADB-PINACA and its metabolites)	10	
Tapentadol	200	
Tramadol	200	<input checked="" type="checkbox"/>
Zolpidem	20	

1. Clinical Laboratory Improvement Amendments (CLIA) 2. Electronic Medical Records (EMR) 3. Health Level Seven International, Inc. (HL7®)

JANT PHARMACAL CORPORATION
LABORATORY & POINT OF CARE DIAGNOSTICS

JANT PHARMACAL CORPORATION
16530 Ventura Blvd, Suite 512
Encino, CA 91436
Tel: 800.676.5565
JantDX.com info@JantDX.com

