

# qTOWER<sup>3</sup> INITIAL PURCHASE PROMOTION



## qTOWER<sup>3</sup> 84

The qTOWER<sup>3</sup> utilizes eight individual light fibers capable of scanning an entire plate in only six seconds and eliminating the need for a passive reference.

High-performance systems for high-throughput applications using real-time PCR in a **384-sample multiwell format**.

State-of-the-art block heating technology and temperature homogeneity across all wells is guaranteed.

**PROMO  
PRICE\***

**\$39,900**

\* Promo price includes the following when RT-PCR panels are added to the analyzer installation:

### EQUIPMENT

qTOWER<sup>3</sup> with all 6 color modules

Free IQ/OQ upon installation

10-year warranty on all optical components

3-year warranty on the remainder of the system

### ★ SERVICE ★

Starting year 2, receive one year of extended Gold Star Service Package for the qTOWER<sup>3</sup>, including Maintenance with OQ: general checkup, inspection & cleaning, validation & performance testing, certification, travel costs, on-site labor, parts and more

### RT-PCR PANELS

Get 20% off on your initial order of pre-plated customized RT-PCR panels

Receive complimentary remote or on-site validation assistance for custom panels

(Up to \$10,000 value)


# INCLUDES 6 COLOR MODULES (\$4,000 value each)



- BLUE (FAM)
- GREEN (JOE, VIC, HEX)
- YELLOW (TAMRA, HEX)
- ORANGE (ROX)
- RED (Cy5)
- NIR1 (Cy5.5)



**GET 20% OFF**  
**PRE-PLATED PANELS**



Customize your target list from the following panels and we will supply you with pre-spotted plates based on your requirements:

## RT-PCR Panels

## Number of Targets Available

Respiratory	.....36
Wound	.....36
Urinary Tract Infection (UTI)	.....36
Antibiotic Resistance	.....24
Nail Fungus	.....14
Sexual Health	.....11

Ready to speak to our technical staff?

**Call (800) 676-5565**