


@accutest[®]

URINE COTININE TEST

CAT#DS-71

40 TESTS/BOX

- **The Accutest[®] Urine Cotinine Test is a rapid, one-step urine drug screen cassette for the simultaneous, qualitative detection of tobacco and nicotine in human urine.**
- **The test detects Cotinine in urine at a cutoff of 200 ng/mL.**



FEATURES:

- One-step drug screen
- Results in 5 minutes
- Easy cassette format - Just add urine
- Disposable
- Built-in control to validate results
- No reagents needed
- Stores at room temperature
- Cut-off level – 200 ng/ml

JANT PHARMACAL CORPORATION
LABORATORY & POINT OF CARE DIAGNOSTICS

accutest is a registered trademark of Jant Pharmacal Corporation.

Jant Pharmacal Corporation
16530 Ventura Blvd., Suite 512
Encino, CA 91436
800.676.5565 Fax: 818.986.0235
info@JantDX.com JantDX.com

CONTACT US FOR MORE INFO

800.676.5565