

accutest®


URINE REAGENT STRIPS

100 TESTS PER BOTTLE

Accutest® Urine Reagent Strips for Urinalysis are in vitro diagnostic test devices that use reagents for qualitative and semi-quantitative urinalysis.

The strips are for single use in professional near patient (point-of-care) facilities and centralized laboratory locations by medical technologists. CLIA-Waived when read visually.

- ✓ HIGH QUALITY TEST STRIPS
- ✓ CLIA-WAIVED WHEN READ VISUALLY
- ✓ FAST QUALITATIVE RESULTS
- ✓ 100 TESTS PER BOTTLE
- ✓ RELIABLE SENSITIVITY
- ✓ EASY TO READ



URS-4 Urine Reagent Strips

CAT.#UA774 - 100 Tests/Bottle
4 Parameter: pH, Glucose, Protein, Ketone.

URS-10 Urine Reagent Strips

CAT.#UA710B - 100 Tests/Bottle
CAT.#UA824B - 24 Bottles/Case
10 Parameter: Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ketone, Bilirubin, Glucose.

URS-11 Urine Reagent Strips

CAT.#UA711 - 100 Tests/Bottle
CAT.#UA924 - 24 Bottles/Case
11 Parameter: Leukocytes, Nitrite, Urobilinogen, Protein, pH, Blood, Specific Gravity, Ascorbic Acid, Ketone, Bilirubin, Glucose.

CONTACT US FOR MORE INFO
800.676.5565

JANT PHARMACAL CORPORATION
LABORATORY & POINT OF CARE DIAGNOSTICS

accutest® is a registered trademark of Jant Pharmacal Corporation.

Jant Pharmacal Corporation
16530 Ventura Blvd., Suite 512
Encino, CA 91436
800.676.5565 Fax: 818.986.0235
info@JantDX.com JantDX.com