

# accutest® ValuPak™ Series

Quality and Convenience at a Lower Cost

## Overview

### Accutest® ValuPak™ iFOBT Test



### Accutest® ValuPak™ Pregnancy Test



### Accutest® ValuPak™ Drug Test Cups

The new ValuPak™ product line offers the same quality and convenience that you have come to expect from the Accutest® brand but at a lower cost per test. The ValuPak™ menu of FDA cleared and CLIA-waived tests will continue to expand in the future, to provide rapid testing solutions for varied clinical practice settings.

## Featured Products



### Accutest® ValuPak™ iFOBT Test

- ✓ Rapid Immunoassay test specific to human hemoglobin
- ✓ Detects hemoglobin concentrations as low as 50 ng/mL
- ✓ High sensitivity, specificity, & accuracy
- ✓ Results within 5 minutes
- ✓ CLIA-Waived

25 Tests per Kit

Cat #	Product Name	CPT Code
CS450	Accutest® ValuPak™ iFOBT Test	81025QW
CS620	Accutest® iFOBT Specimen Collection Kits	
CS430	Accutest® iFOBT Collection Tubes	

### Accutest® ValuPak™ Pregnancy Test

- ✓ For the qualitative detection of hCG in urine specimens
- ✓ Detects hCG concentrations as low as 25 mIU/mL
- ✓ High sensitivity, specificity, & accuracy
- ✓ One-step test, no reagents needed
- ✓ Easy to use cassette format
- ✓ Results within 5 minutes
- ✓ CLIA-Waived

40 Tests per Kit

Cat #	Product Name	CPT Code
PF453	Accutest® ValuPak™ Pregnancy Test	82274QW

### Accutest® ValuPak™ Drug Test Cups

- ✓ Screens for the presence of multiple drugs
- ✓ High sensitivity, specificity, & accuracy
- ✓ Available in 6 and 12 panel formats
- ✓ Adulteration testing available
- ✓ Results within 5 minutes
- ✓ CLIA-Waived

25 Tests per Kit

Cat #	Product Name	Configuration	CPT Code
DS412	Accutest® ValuPak™ Drug Test Cup 12	AMP/BAR/BUP/BZO/COC/MDMA/MET/MOP/MTD/OXY/PCP/THC	80305QW
DS432	Accutest® ValuPak™ Drug Test Cup 12 with Adulterants	AMP/BAR/BUP/BZO/COC/MDMA/MET/MOP/MTD/OXY/PCP/THC + (OXI/SG/PH)	80305QW
DS487	Accutest® ValuPak™ Drug Test Cup 6	AMP/COC/MET/OPI/PCP/THC	80305QW