



H. PYLORI UREASE TEST CASSETTE



**For Rapid Detection of
Helicobacter Pylori
in Gastric Biopsies**

ACCUTEST H. PYLORI UREASE TEST FEATURES

NO REFRIGERATION REQUIRED

Can be stored at room temperature.

SIMPLE & CONVENIENT

May be used in any CLIA-waived endoscopy unit with minimal training. Simply place biopsy in the test well, close and wait for results.

RAPID RESULTS

Positive results in as little as 1 minute, negative results in 1 hour.

ACCURATE

>99% Sensitivity and Specificity.*

DRY TEST

No incubation or reagents required.

COST EFFECTIVE

Low cost per test.

3 YEAR SHELF LIFE

Good stability, superior shelf life.

CPT CODE

87077QW

*Refer to the Instructions for Use (IFU) for detailed information and studies.



CATALOG

ID750

IDC73

DESCRIPTION

Accutest H. pylori Urease Test

Accutest Urease/Glycerin Controls Solution

UNITS

50 tests/box (control included)

2 mL dropper bottle

**JANT PHARMACAL
CORPORATION**
LABORATORY & POINT OF CARE DIAGNOSTICS

accutest® is a registered trademark of Jant Pharmacal Corporation.

800.676.5565

info@jantdx.com

jantdx.com